XP Screens Crystal Screening with the Protein Crystal Glue TEW





The **XP Screen** is a convenient initial screen that has successfully induced protein crystallization with low TEW concentrations such as 1 mM^[5,6]. In some cases however, higher concentration of 5 or 10 mM TEW are needed.

The **XP Up Screen** is an upgrade of the well-established XP Screen. It contains 96 of the most prominent screening solutions that are long-term stable in the presence of up to 10 mM TEW.

Advantages of XP Screens

- Convenient initial screening with up to 10 mM TEW
- Improved crystal quality and resolution
- Tight crystal contacts
- New crystal forms (e.g to provide access to binding pockets)
- Direct phasing



Figure from [1].

Product	CatNo.	Amount
XP Screen	CS-350	96 screening conditions with 1 mM TEW
XP Up Screen	CS-351	96 screening conditions to be upgraded with up to 10 mM TEW
Anderson-Evans polyoxotungstate (TEW)	X-TEW-5	Additive for optimization

Crystallography & Cryo-EM

XP Screens Crystal Screening with the Protein Crystal Glue TEW





The Anderson-Evans Polyoxotungstate [TeW₆O₂₄]⁶⁻ (TEW) is a universal and flexible crystallization additive^[1] that is integrated in our XP Screens. It was shown to improve crystal quality and resolution by:

- binding to the protein surface and forming tight crystal lattice contacts^[3-7]
- acting as linker in various orientations and thereby creating either smaller (PDB: 40UA) or larger (PDB: 4PHI) protein-protein distances^[3,4]
- heterogeneous crystallization of two different protein forms in one single crystal (PDB: 4OUA)^[4]
- inducing new crystal forms with access to the active site (PDB: 6ZDK, 6ZDL, 6ZFA)^[5]
- covalently binding and structurally adapting to fit into protein molecules (PDB: 4Z12, 4Z13)^[2]



Figure 1: Covalent binding of Tungsten to carboxylic oxygen atoms of glutamic acid (PDB: 4Z12, 4Z13).^[2] Figure from [1].



TEW Properties

- Centrosymmetric, inorganic cluster
- Disk-shaped ion [TeW₆O₂₄]⁶⁻
- Dimensions: 9 x 9 x 3 Å³
- MW: 1.615 g/mol (ion)
- Highly soluble in aqueous solutions (100 mM) and stable over a wide pH range
- 6 tungsten atoms provide strong anomalous scattering signal for phasing

Figure 2: Chemical structure of TEW. Figure from [7].

References:

- [1] Bijelic et al. (2017) Ten Good Reasons for the Use of the Tellurium-Centered Anderson-Evans Polyoxotungstate in Protein Crystallography. Acc. Chem. Res. 50:1441.
- [2] Molitor et al. (2016) In situ formation of the first proteinogenically functionalized [TeW602402(Glu)]7- structure reveals unprecedented chemical and geometrical features of the Anderson-type cluster. Chem. Commun. 52:12286.
- [3] Bijelic *et al.* (2015) Hen Egg-White Lysozyme Crystallisation: Protein Stacking and Structure Stability Enhanced by a Tellurium(VI)-Centred Polyoxotungstate. *ChemBioChem* **16**:233.
- [4] Mauracher et al. (2014) Latent and active abPPO4 mushroom tyrosinase cocrystallized with hexatungstotellurate(VI) in a single crystal. Acta Cryst. D 70:2301.
- [5] Sobala et al. (2020) Structure of human endo-α-1,2-mannosidase (MANEA), an antiviral host-glycosylation target. PNAS 117 (47):29595.
- [6] Ames et al. (2020) Identifying a Molecular Mechanism That Imparts Species-Specific Toxicity to YoeB Toxins. Front. Microbiol. 11:959.
- [7] Mac Sweeney et al. (2018) The crystallization additive hexatungstotellurate promotes the crystallization of the HSP70 nucleotide binding domain into two different crystal forms. PLOS one **13 (6)**: e0199639.



Jena Bioscience GmbH Loebstedter Str. 71 07749 Jena Germany Phone +49(0)3641-6285 000 Fax +49(0)3641-6285 100 info@jenabioscience.com www.jenabioscience.com



IFTA AG Certified QMS and EMS according to DIN EN ISO 9001 and DIN EN ISO 14001 Reg.-No.: ICV03597 034 and ICV03597 534

Crystallography & Cryo-EM