

**LEXSY host P10**

glycerol stocks
for constitutive expression

Cat. No.	Amount
LT-101	3 x 1,6 ml

For research use only! Not intended for human or animal diagnostic or therapeutic uses.

Shipping: Shipped on dry ice

Storage Conditions: Store at -80 °C

Additional Storage Conditions: Avoid freeze/thaw cycles

Upon arrival the glycerol stocks must be stored at -80 °C or inoculated into LEXSY BHI.

Shelf Life: 12 months

Description:

The *Leishmania tarentolae* laboratory strain P10 is used as host strain for the constitutive LEXSY expression vectors.

Content:

3 vials with 1.6 ml each of frozen glycerol stocks of LEXSY host P10 (*Leishmania tarentolae* laboratory strain P10).

These stocks can be stored at -80 °C for at least 1 year. For reactivation see below.

Organism:

Leishmania tarentolae, laboratory strain P10

Biosafety level:

1, Non-pathogenic for mammals

Source:

Tarentola annularis

Reactivation of LEXSY host:

Thaw glycerol stock on ice and inoculate the entire content of the vial into 10 ml of LEXSY BHI medium (Cat.-No. ML-411). Incubate at 26 °C and dilute as required.

Preparation of LEXSY BHI growth medium:

Dissolve 37 g/l LEXSY BHI powder (Cat.-No. ML-412) in deionized water and autoclave for 15 min at 121 °C. Add Hemin and PenStrep. Store at 4 °C and use within two weeks.

Preparation of glycerol stocks:

Add 1.2 ml of growing culture (ca. 6×10^7 cells/ml) to one vial with 0.4 ml 80 % glycerol, mix, incubate 10 min at RT, 1 h on ice and overnight at -20 °C. Transfer to -80 °C. Strains can be stored this way for several years.

Selected References:

[1] Zauner *et al.* (2018) Structural Analyses of *Arabidopsis thaliana* Legumain γ Reveal Differential Recognition and Processing of Proteolysis and Ligation Substrates. *J Biol Chem* **293**: 8934