

JBScreen PACT++ 4

Cat.-No: CS-164

Number	Precipitant	Buffer	pH	Salt
A 1	20% w/v PEG 3350	100 mM Bis-Tris Propane	7.5	200 mM Sodium Fluoride
A 2	20% w/v PEG 3350	100 mM Bis-Tris Propane	7.5	200 mM Sodium Bromide
A 3	20% w/v PEG 3350	100 mM Bis-Tris Propane	7.5	200 mM Sodium Iodide
A 4	20% w/v PEG 3350	100 mM Bis-Tris Propane	7.5	200 mM Potassium Thiocyanate
A 5	20% w/v PEG 3350	100 mM Bis-Tris Propane	7.5	200 mM Sodium Nitrate
A 6	20% w/v PEG 3350	100 mM Bis-Tris Propane	7.5	200 mM Sodium Formate
B 1	20% w/v PEG 3350	100 mM Bis-Tris Propane	7.5	200 mM Sodium Acetate
B 2	20% w/v PEG 3350	100 mM Bis-Tris Propane	7.5	200 mM Sodium Sulfate
B 3	20% w/v PEG 3350	100 mM Bis-Tris Propane	7.5	200 mM Potassium/Sodium Tartrate
B 4	20% w/v PEG 3350	100 mM Bis-Tris Propane	7.5	200 mM Sodium/Potassium Phosphate
B 5	20% w/v PEG 3350	100 mM Bis-Tris Propane	7.5	200 mM Sodium Citrate
B 6	20% w/v PEG 3350	100 mM Bis-Tris Propane	7.5	200 mM Sodium Malonate
C 1	20% w/v PEG 3350	100 mM Bis-Tris Propane	8.5	200 mM Sodium Fluoride
C 2	20% w/v PEG 3350	100 mM Bis-Tris Propane	8.5	200 mM Sodium Bromide
C 3	20% w/v PEG 3350	100 mM Bis-Tris Propane	8.5	200 mM Sodium Iodide
C 4	20% w/v PEG 3350	100 mM Bis-Tris Propane	8.5	200 mM Potassium Thiocyanate
C 5	20% w/v PEG 3350	100 mM Bis-Tris Propane	8.5	200 mM Sodium Nitrate
C 6	20% w/v PEG 3350	100 mM Bis-Tris Propane	8.5	200 mM Sodium Formate
D 1	20% w/v PEG 3350	100 mM Bis-Tris Propane	8.5	200 mM Sodium Acetate
D 2	20% w/v PEG 3350	100 mM Bis-Tris Propane	8.5	200 mM Sodium Sulfate
D 3	20% w/v PEG 3350	100 mM Bis-Tris Propane	8.5	200 mM Potassium/Sodium Tartrate
D 4	20% w/v PEG 3350	100 mM Bis-Tris Propane	8.5	200 mM Sodium/Potassium Phosphate
D 5	20% w/v PEG 3350	100 mM Bis-Tris Propane	8.5	200 mM Sodium Citrate
D 6	20% w/v PEG 3350	100 mM Bis-Tris Propane	8.5	200 mM Sodium Malonate