

Pentaerythritol ethoxylate (3/4 EO/OH)

Pentaerythritol ethoxylate (3/4 EO/OH), 60 % w/v

PEE 270

Cat. No.	Amount
CSS-372	100 ml

Application: Reagent for optimization of crystallization conditions for proteins.

JBScreen Single Stocks - PEE 270 is a ready-made solution.

The single stock is provided as a 60% w/v sterile filtered stock solution and supplied in 100 ml aliquots.

Molecular Formula: $C[CH_2(OCH_2CH_2)_nOH]_4$

CAS Number: 30599-15-6

EC Number: 500-071-2

Storage: 4°C

Expiry: 12 months