

Pentaerythritol propoxylate (17/8 PO/OH)

Pentaerythritol propoxylate (17/8 PO/OH), 60 % w/v

PEP 629

Cat. No.	Amount
CSS-376	100 ml

Application: Reagent for optimization of crystallization conditions for proteins.

JBScreen Single Stocks - PEP 629 is a ready-made solution.

The single stock is provided as a 60% w/v sterile filtered stock solution and supplied in 100 ml aliquots.

Molecular Formula: $C[CH_2[OCH_2CH(CH_3)]_nOH]_4$

CAS Number: 9051-49-4

EC Number: 500-030-9

Storage: 4°C

Expiry: 12 months