

## Pentaerythritol propoxylate (17/8 PO/OH)

Pentaerythritol propoxylate (17/8 PO/OH), 100 % v/v

PEP 629

Cat. No.	Amount
CSS-377	100 ml

**Application:** Reagent for optimization of crystallization conditions for proteins.

**JBScreen Single Stocks - PEP 629** is a ready-made solution.

The single stock is provided as a 100% v/v sterile filtered stock solution and supplied in 100 ml aliquots.

**Molecular Formula:**  $C[CH_2[OCH_2CH(CH_3)]_nOH]_4$

**CAS Number:** 9051-49-4

**EC Number:** 500-030-9

**Storage:** 4°C

**Expiry:** 12 months