



## MicroTools™ Kit 1

20 soft, flexible tools to simplify your sample handling tasks




**MicroSpoons™**

Push, detach from surfaces, separate from a cluster, lift, transfer. Push down to flatten tip and slide under sample.




**MicroRuler™**

Scale in microns.




**MicroChisels™**

Push, detach from surfaces, separate from a cluster, lift, transfer. Push down to flatten tip and slide under sample.




**MicroSieves™**

Sieve, lift, transfer and soak.




**MicroTip™**

Poke, separate and lift.




**MicroLassos™**

Lasso and pull to detach from surfaces.



**MicroSaws™**

Cut and remove protein skins, gels, lipids.



**MicroGrippers™**

Push down onto your sample to firmly capture it for room-temperature crystallography.