


MicroTools™ Kit 3


20 soft, flexible tools for linear, radial and angular measurements
Available in both standard and X-ray fluorescent versions



1-5

Horizontal MicroRuler™


500 μm wide with 25 μm increments
For linear measurements



6-10

Vertical MicroRuler™


1000 μm wide with 25 μm increments
For linear measurements



11-15

Large MicroProtractor™

1000 μm wide with 50 μm and 15° increments
For radial and angular measurements



16-20

Small MicroProtractor™

320 μm wide with 20 μm and 45° increments
For radial and angular measurements