

## Products by Feature:

	MicroMounts	MicroLoops	MicroLoops HT	MicroLoops LD	MicroLoops E	MicroMeshes	MicroGrippers
<i>Aperture Sizes (μm)</i>	10-600	20-300	20-300	20-300	150-700 Elongated	10-50 Grid	50-300 Thin Fingers
<i>Thickness (μm)</i>	10-12.5	10-12.5	18	Dual Thickness	10-12.5	10-12.5	10-12.5
<i>Background Scatter</i>	✓✓	✓✓	✓	✓✓✓	✓✓✓	✓✓	✓✓
<i>Rigidity</i>	✓✓	✓✓	✓✓✓	✓✓	✓✓	✓✓	✓✓
<i>Durability</i>	✓✓	✓✓	✓✓✓	✓✓	✓✓	✓✓	✓✓
<i>Sample Visibility</i>	✓✓	✓✓	✓✓	✓✓✓	✓✓✓	✓	✓✓
<i>Ease of Liquid Removal</i>	✓✓✓	✓✓	✓✓	✓✓	✓✓	✓	✓✓

✓✓✓	Excellent choice for required features and applications
✓✓	Good choice for required features and applications
✓	Adequate choice for required features and applications

## Products by Application:

	MicroMounts	MicroLoops	MicroLoops HT	MicroLoops LD	MicroLoops E	MicroMeshes	MicroGrippers
<i>Protein</i>	✓✓	✓✓	✓✓	✓✓✓	✓✓✓	✓✓✓	✓✓
<i>Small Molecule</i>	✓✓	✓✓	✓✓	✓✓	✓✓	✓✓	✓✓
<i>Microcrystal</i>	✓✓	✓✓	✓	✓✓✓	✓✓✓	✓✓✓	✓
<i>Thin Plates</i>	✓✓	✓✓	✓	✓✓	✓✓	✓✓✓	✓✓✓
<i>Needles and Rods</i>	✓✓	✓✓	✓	✓✓	✓✓✓	✓✓	✓✓✓
<i>Room Temperature</i>	✓✓	✓✓	✓✓	✓✓	✓✓	✓✓	✓✓✓

✓✓✓	Excellent choice for required features and applications
✓✓	Good choice for required features and applications
✓	Adequate choice for required features and applications