


2'-I-ADP

2'-Iodo-Adenosine-5'-diphosphate, Sodium salt

Cat. No.	Amount
NU-109S	30 Units
NU-109L	150 Units



Cat. No.: NU-109

Molecular Formula: C₁₀H₁₄N₅O₉P₂I (free acid)

Molecular Weight: 537.10 (free acid)

Purity: > 95%, clear aqueous solution, pH 7.5

Storage conditions:

Short term exposure (up to 1 week cumulative) to ambient temperature possible. Long term storage at < -20°C. If stored as recommended, Jena Bioscience guarantees optimal performance of this product for 12 months after date of delivery.

For research use only!

1 unit = 1 µl of a 10 mM solution