


## 2'Br-ADP

2'-Bromo-Adenosine-5'-diphosphate, Sodium salt

Cat. No.	Amount
NU-110S	30 Units
NU-110L	150 Units



**Cat. No.:** NU-110

**Molecular Formula:**  $C_{10}H_{14}N_5O_9P_2Br$  (free acid)

**Molecular Weight:** 490.10 (free acid)

**Purity:** > 95%, clear aqueous solution, pH 7.5

**Storage conditions:**

Short term exposure (up to 1 week cumulative) to ambient temperature possible. Long term storage at < -20°C. If stored as recommended, Jena Bioscience guarantees optimal performance of this product for 12 months after date of delivery.

**For research use only!**

1 unit = 1  $\mu$ l of a 10 mM solution