


8I-ADP

8-Iodo-Adenosine-5'-diphosphate, Sodium salt

Cat. No.	Amount
NU-111S	50 Units
NU-111L	250 Units



Cat. No.: NU-111

Molecular Formula: $C_{10}H_{14}N_5O_{10}P_2I$ (free acid)

Molecular Weight: 553.10 (free acid)

Purity: > 95%, clear aqueous solution, pH 7.5

Storage conditions:

Short term exposure (up to 1 week cumulative) to ambient temperature possible. Long term storage at $< -20^{\circ}C$. If stored as recommended, Jena Bioscience guarantees optimal performance of this product for 12 months after date of delivery.

For research use only!

1 unit = 1 μ l of a 10 mM solution