


8I-AppNHp (8I-AMPPNP)

8-Iodo-adenosine-5'-[(β,γ)-imido]triphosphate, Sodium salt

Cat. No.	Amount
NU-115S	20 Units
NU-115L	100 Units



Cat. No.: NU-115

Molecular Formula: C₁₀H₁₆N₆O₁₂P₃I (free acid)

Molecular Weight: 632.09 (free acid)

Purity: > 95%, clear aqueous solution, pH 7.5

Storage conditions:

Short term exposure (up to 1 week cumulative) to ambient temperature possible. Long term storage at < -20°C. If stored as recommended, Jena Bioscience guarantees optimal performance of this product for 6 months after date of delivery.

For research use only!

1 unit = 1 μ l of a 10 mM solution