


## 8I-GTP

### 8-Iodo-guanosine-5'-triphosphate, Sodium salt

Cat. No.	Amount
NU-125S	30 Units
NU-125L	150 Units



**Cat. No.:** NU-125

**Molecular Formula:**  $C_{10}H_{15}N_5O_{14}P_3I$  (free acid)

**Molecular Weight:** 649.07 (free acid)

**Purity:** > 95%, clear aqueous solution, pH 7.5

**Storage conditions:**

Short term exposure (up to 1 week cumulative) to ambient temperature possible. Long term storage at < -20°C. If stored as recommended, Jena Bioscience guarantees optimal performance of this product for 12 months after date of delivery.

**For research use only!**

1 unit = 1  $\mu$ l of a 10 mM solution