


8Br-GMP

8-Bromo-Guanosine-5'-monophosphate, Sodium salt

Cat. No.	Amount
NU-144S	100 Units
NU-144L	500 Units



Cat. No.: NU-144

Molecular Formula: $C_{10}H_{13}N_5O_8PBr$ (free acid)

Molecular Weight: 442.11 (free acid)

Purity: > 95%, clear aqueous solution, pH 7.5

Storage conditions:

Short term exposure (up to 1 week cumulative) to ambient temperature possible. Long term storage at $< -20^{\circ}C$. If stored as recommended, Jena Bioscience guarantees optimal performance of this product for 12 months after date of delivery.

* 1 unit = 1 μ l of a 10 mM solution

For research use only!