


## 2'-I-AMP

### 2'-Iodo-Adenosine-5'-monophosphate, Sodium salt

Cat. No.	Amount
NU-147S	50 Units
NU-147L	250 Units



**Cat. No.:** NU-147

**Molecular Formula:** C<sub>10</sub>H<sub>13</sub>N<sub>5</sub>O<sub>6</sub>PI (free acid)

**Molecular Weight:** 457.12 (free acid)

**Purity:** > 95%, clear aqueous solution, pH 7.5

**Storage conditions:**

Short term exposure (up to 1 week cumulative) to ambient temperature possible. Long term storage at < -20°C. If stored as recommended, Jena Bioscience guarantees optimal performance of this product for 12 months after date of delivery.

\* 1 unit = 1 µl of a 10 mM solution

**For research use only!**