


2'-Cl-ATP

2'-Chloro-Adenosine-5'-triphosphate, Sodium salt

Cat. No.	Amount
NU-148S	50 Units
NU-148L	250 Units



Cat. No.: NU-148

Molecular Formula: $C_{10}H_{15}N_5O_{12}P_3Cl$ (free acid)

Molecular Weight: 525.62 (free acid)

Purity: > 95%, clear aqueous solution, pH 7.5

Storage conditions: Short term exposure (up to 1 week cumulative) to ambient temperature possible. Long term storage at $< -20^{\circ}C$. If stored as recommended, Jena Bioscience guarantees optimal performance of this product for 12 months after date of delivery.

1 unit = 1 μ l of a 10 mM solution

For research use only!