




pLEXY_I-ble3

Integrative inducible expression vector
Antibiotic selection of transfectants with LEXSY Bleo

Cat. No.	Amount
EGE-244	5 µg

Description:

Inducible LEXSY expression vector pLEXY_I-ble3 with the *ble* marker gene (*bleomycin resistance gene* from *S. hindustans*), allowing selection of recombinant LEXSY strains with the antibiotic LEXSY Bleo (Cat.-No. AB-103). Designed for inducible expression of target genes in LEXSY host T7-TR (Cat.-No. LT-110).



For research use only! Not intended for human or animal diagnostic or therapeutic uses.

Shipping: shipped on gel packs

Storage Conditions: Store at -20 °C

Shelf Life: 12 months