




## pLEXY-ble2.1

Integrative constitutive expression vector  
Antibiotic selection of transfectants with LEXSY Bleo

Cat. No.	Amount
EGE-271	5 µg

### Description:

General expression vector pLEXY-ble2.1 with the *bleomycin resistance gene* from *S. hindustans* allowing selection of recombinant LEXSY strains with the antibiotic LEXSY Bleo (Cat.-No. AB-103).



**For research use only!** Not intended for human or animal diagnostic or therapeutic uses.

**Shipping:** shipped on gel packs

**Storage Conditions:** Store at -20 °C

**Shelf Life:** 12 months