




pLEXY-hyg2.1

integrative constitutive expression vector
antibiotic selection of transfectants with LEXSY Hygro

Cat. No.	Amount
EGE-272	5 µg

Description:

General expression vector pLEXY-hyg2.1 with the *hyg* marker gene (encoding Hygromycin phosphotransferase, HPT) allowing selection of recombinant LEXSY strains with the antibiotic LEXSY Hygro (Cat.-No. AB-104).



For research use only! Not intended for human or animal diagnostic or therapeutic uses.

Shipping: shipped on gel packs

Storage Conditions: Store at -20 °C

Shelf Life: 12 months