




## pLEXY-neo2.1

Integrative constitutive expression vector  
Antibiotic selection of transfectants with LEXSY Neo

Cat. No.	Amount
EGE-273	5 µg

### Description:

General expression vector pLEXY-neo2.1 with *neo* marker gene encoding aminoglycoside phospho-transferase 3' (I) (APH(3')I) allowing selection of recombinant LEXSY strains with the antibiotic LEXSY Neo (Cat.-No. AB-105).



**For research use only!** Not intended for human or animal diagnostic or therapeutic uses.

**Shipping:** shipped on gel packs

**Storage Conditions:** Store at -20 °C

**Shelf Life:** 12 months