



## pLEXY-sat2.1


Integrative constitutive expression vector

Antibiotic selection of transfectants with Nourseothricin (NTC, clonNAT)

Cat. No.	Amount
EGE-274	5 µg

### Description:

General expression vector pLEXY-sat2.1 with *sat* marker gene (encoding streptothricine acetyltransferase) allowing selection of recombinant LEXSY strains with the antibiotic LEXSY NTC (Cat.-No. AB-101 and Cat.-No. AB-102).



**For research use only!** Not intended for human or animal diagnostic or therapeutic uses.

**Shipping:** shipped on gel packs

**Storage Conditions:** Store at -20 °C

**Shelf Life:** 12 months