




## pLEXY-pac2.1

Integrative constitutive expression vector  
Antibiotic selection of transfectants with LEXSY Puro

Cat. No.	Amount
EGE-276	5 µg

### Description:

General expression vector pLEXY-pac2.1 with *pac* marker gene (puromycin N-acetyltransferase from *Streptomyces alboniger*) allowing selection of recombinant LEXSY strains with the antibiotic LEXSY Puro (Cat.-No. AB-108).



**For research use only!** Not intended for human or animal diagnostic or therapeutic uses.

**Shipping:** shipped on gel packs

**Storage Conditions:** Store at -20 °C

**Shelf Life:** 12 months