

**LEXSY in vitro Translation Cell Extract**

for cell-free protein synthesis

Cat. No.	Amount
EGE-260	60 reactions (2 x 500 µl)

**For research use only!** Not intended for human or animal diagnostic or therapeutic uses.

**Shipping:** shipped on dry ice

**Storage Conditions:** store at -80 °C

**Shelf Life:** 12 months

**Description:**

The LEXSY cell extract for *in vitro* translation contains functional ribosomes and all essential components of the eukaryotic translation and folding machinery. For target mRNA generation the heterologous T7 RNA polymerase was added to the extracts.

**Please note:**

EGE-260 is part of the LEXSY *in vitro* Translation Kit available as plasmid based (EGE-2002-15) or PCR based (EGE-2010-15) cell-free protein synthesis kit.