

**Pen-Strep**

sterile 200x stock solution of penicillin and streptomycin

Cat. No.	Amount
ML-105S	5 ml (for 1 L medium)
ML-105L	25 ml (for 5 L medium)
ML-105XL	50 ml (for 10 L medium)
ML-105XXL	250 ml (for 50 L medium)

For research use only! Not intended for human or animal diagnostic or therapeutic uses.

Shipping: shipped on gel packs

Storage Conditions: store at -20 °C

Shelf Life: 12 months

Description:

Contains 10.000 units of penicillin (base) and 10.000 µg of streptomycin (base)/ml as penicillin G sodium salt and streptomycin sulfate in 0.85% saline (200x stock solution).

Preparation of *Leishmania* cultivation media:

Add to 1 l medium:

5 ml 200x Pen-Strep stock solution.