

Wizard Classic 2 tube kit

M-R-1009531

Tube	Precipitation Reagent	Buffer	Salt
1	10% (w/v) PEG 3000	100 mM Sodium acetate/ Acetic acid pH 4.5	200 mM Zinc acetate
2	35% (v/v) MPD	100 mM MES/ Sodium hydroxide pH 6.0	200 mM Lithium sulfate
3	20% (w/v) PEG 8000	100 mM Tris base/ Hydrochloric acid pH 8.5	200 mM Magnesium chloride
4	2000 mM Ammonium sulfate	100 mM Sodium cacodylate/ Hydrochloric acid pH 6.5	200 mM Sodium chloride
5	20% (v/v) 1,4-butanediol	100 mM HEPES/ Sodium hydroxide pH 7.5	200 mM Sodium chloride
6	10% (v/v) 2-propanol	100 mM Sodium phosphate dibasic/ Citric acid pH 4.2	200 mM Lithium sulfate
7	30% (w/v) PEG 3000	100 mM Tris base/ Hydrochloric acid pH 7.0	200 mM Sodium chloride
8	10% (w/v) PEG 8000	100 mM Potassium phosphate monobasic/ Sodium phosphate dibasic pH 6.2	200 mM Sodium chloride
9	2000 mM Ammonium sulfate	100 mM Sodium phosphate dibasic/ Citric acid pH 4.2	
10	1000 mM Ammonium phosphate dibasic	100 mM Tris base/ Hydrochloric acid pH 8.5	
11	10% (v/v) 2-propanol	100 mM Sodium cacodylate/ Hydrochloric acid pH 6.5	200 mM Zinc acetate
12	30% (v/v) PEG 400	100 mM Sodium cacodylate/ Hydrochloric acid pH 6.5	200 mM Lithium sulfate
13	15% (v/v) Reagent alcohol	100 mM Sodium citrate/ Citric acid pH 5.5	200 mM Lithium sulfate
14	20% (w/v) PEG 1000	100 mM Potassium phosphate monobasic/ Sodium phosphate dibasic pH 6.2	200 mM Sodium chloride
15	1260 mM Ammonium sulfate	100 mM HEPES/ Sodium hydroxide pH 7.5	
16	1000 mM Sodium citrate tribasic	100 mM CHES/ Sodium hydroxide pH 9.5	
17	2500 mM Sodium chloride	100 mM Tris base/ Hydrochloric acid pH 7.0	200 mM Magnesium chloride
18	20% (w/v) PEG 3000	100 mM Tris base/ Hydrochloric acid pH 7.0	200 mM Calcium acetate
19	1600 mM Sodium phosphate monobasic/ 400 mM Potassium phosphate dibasic	100 mM Sodium phosphate dibasic/ Citric acid pH 4.2	
20	15% (v/v) Reagent alcohol	100 mM MES/ Sodium hydroxide pH 6.0	200 mM Zinc acetate
21	35% (v/v) MPD	100 mM Sodium acetate/ Acetic acid pH 4.5	
22	10% (v/v) 2-propanol	100 mM Imidazole/ Hydrochloric acid pH 8.0	
23	15% (v/v) Reagent alcohol	100 mM HEPES/ Sodium hydroxide pH 7.5	200 mM Magnesium chloride
24	30% (w/v) PEG 8000	100 mM Imidazole/ Hydrochloric acid pH 8.0	200 mM Sodium chloride
25	35% (v/v) MPD	100 mM HEPES/ Sodium hydroxide pH 7.5	200 mM Sodium chloride
26	30% (v/v) PEG 400	100 mM CHES/ Sodium hydroxide pH 9.5	
27	10% (w/v) PEG 3000	100 mM Sodium cacodylate/ Hydrochloric acid pH 6.5	200 mM Magnesium chloride
28	20% (w/v) PEG 8000	100 mM MES/ Sodium hydroxide pH 6.0	200 mM Calcium acetate
29	1260 mM Ammonium sulfate	100 mM CHES/ Sodium hydroxide pH 9.5	200 mM Sodium chloride
30	20% (v/v) 1,4-butanediol	100 mM Imidazole/ Hydrochloric acid pH 8.0	200 mM Zinc acetate
31	1000 mM Sodium citrate tribasic	100 mM Tris base/ Hydrochloric acid pH 7.0	200 mM Sodium chloride
32	20% (w/v) PEG 1000	100 mM Tris base/ Hydrochloric acid pH 8.5	
33	1000 mM Ammonium phosphate dibasic	100 mM Sodium citrate/ Citric acid pH 5.5	200 mM Sodium chloride
34	10% (w/v) PEG 8000	100 mM Imidazole/ Hydrochloric acid pH 8.0	
35	800 mM Sodium phosphate monobasic/ 1200 mM Potassium phosphate dibasic	100 mM Sodium acetate/ Acetic acid pH 4.5	
36	10% (w/v) PEG 3000	100 mM Sodium phosphate dibasic/ Citric acid pH 4.2	200 mM Sodium chloride
37	1000 mM Potassium sodium tartrate	100 mM Tris base/ Hydrochloric acid pH 7.0	200 mM Lithium sulfate
38	2500 mM Sodium chloride	100 mM Sodium acetate/ Acetic acid pH 4.5	200 mM Lithium sulfate
39	20% (w/v) PEG 8000	100 mM CAPS/ Sodium hydroxide pH 10.5	200 mM Sodium chloride
40	20% (w/v) PEG 3000	100 mM Imidazole/ Hydrochloric acid pH 8.0	200 mM Zinc acetate
41	2000 mM Ammonium sulfate	100 mM Tris base/ Hydrochloric acid pH 7.0	200 mM Lithium sulfate
42	30% (v/v) PEG 400	100 mM HEPES/ Sodium hydroxide pH 7.5	200 mM Sodium chloride
43	10% (w/v) PEG 8000	100 mM Tris base/ Hydrochloric acid pH 7.0	200 mM Magnesium chloride
44	20% (w/v) PEG 1000	100 mM Sodium cacodylate/ Hydrochloric acid pH 6.5	200 mM Magnesium chloride
45	1260 mM Ammonium sulfate	100 mM MES/ Sodium hydroxide pH 6.0	
46	1000 mM Ammonium phosphate dibasic	100 mM Imidazole/ Hydrochloric acid pH 8.0	200 mM Sodium chloride
47	2500 mM Sodium chloride	100 mM Imidazole/ Hydrochloric acid pH 8.0	200 mM Zinc acetate
48	1000 mM Potassium sodium tartrate	100 mM MES/ Sodium hydroxide pH 6.0	